




Aspen Quality Sacral Stabilization

The new Evergreen™ SI delivers the quality of design and durability that you've come to expect from Aspen.

Code L0621 Approved



Low Profile and Lightweight

Durable two-ply construction with breathable mesh fabric makes the Evergreen SI comfortable for extended wear.

Bi-Lateral Pull

Two easy-to-locate pull tabs provide maximum compression without the spinal rotation that can occur with a single pull tab system.

Effective Compression

Utilizes Aspen's proven SlickTrack™ mechanically advantaged tightening system. Easy to pull regardless of patient strength.



Non-Slip Strips

Silicone strips on the inside of the brace hold fast and prevent migration especially while moving and sitting.



Product Information

| | |
|----------------------|---------|
| Evergreen SI, Small | 36-2102 |
| Evergreen SI, Medium | 36-2104 |
| Evergreen SI, Large | 36-2106 |



Life Changing Spinal Innovation



The Aspen Advantage Symbol

An innovative feature available only from Aspen